Non-image forming effects of light


Lecture in chronobiological phototherapy for doctors. Prague 5


Dipl. Ing. Antonín Fuksa


February 10, 2016

















A novel photoreceptor (i)

- ipRGCs
 <u>Intrinsically photosensitive Retinal Ganglion Cells</u>
- 1923 (Keeler)
- 1991 (Foster, Provencio)
- 1998 (Provencio)
- 2002 (Hattar, Berson)
- 2007 (Zaidi, Zeisler, Foster)


A novel photoreceptor (ii)

Non-visual (non-image forming) function:


to signalize light level to the brain.

- Synchronization of sleep/wake cycle.
 Circadian phase shift. Long term.
- Activation of the organism by light.
 Melatonin Supression. Short term.


Non-image forming effects of light




Spectral response of ipRGCs



ipRGCs distribution on retina




Phase response curve




Zeitzer, JPhys 2000 [1]

Annual cycle of light (50th parallel)



Night length and melatonin



Research in light therapy

- History (classics -> Rosenthal).
- Elderly care (Zumtobel).
- Patient room lighting (Philips, Xal).
- VFN Prague (NASLI).

- Chronobiological treatment of depression.
- Additional treatment for dementia.

Chronobiological phototherapy

- Timing of application.
- Intensity at eyes × length of exposure.
- Contents of the activating blue component.

- Acceptability of light.
- Light dynamics during the day and year.
- Individual factors.

Case study - NASLI SunSun

- Desk luminaires
 10 000 lx @ 0,2-0,3 m
- For personal use SAD and general depression

- $T_c = 6500 \text{ K}$
- P = 72 and 110 W



Case study - NASLI Evidens

- Pendant luminaires
- Patient+lunch room
- E = 30 3000 Ix
- $T_{cp} = 6500 \text{ K}$

VFN Prague

2012





Case study - NASLI ADS2max

- Mobile lumin. to patient's bed
- h = 2,35 m (1,8 m transp.)
- E = 0-2000 Ix
- $T_c = 1800-6500 \text{ K}$
- Autonomous and research modes
- VFNPrague





Autonomous mode of ADS2max


Simple control with remote controller.

• Programs:


Morning (+evening) bright light therapy. All-day phototherapy ("light enhancement"). Down/dusk simulation + combinations. Aid for Wake Therapy.

Autonomous programs

All-day phototherapy




Bright light therapy



Research mode of ADS2max

Maximum flexibility

- Doctor uploads a tailored light schedule on an SD card.
- The luminaire automatically controls the intensity and tone of light depending on the schedule and precision built in clock.



Statutory health insurance

- Chronobiological treatment including phototerapy is a service no.35115 according to Public notice no. 326 / 2014 Coll of Laws.
- In a specialized facility only.
- During patient hospitalization only.
- Competent staff.
- Suitable technical equipment.

Our references

- Psychiatrical clinic 1. LF UK a VFN, Prague
- Alzheimercentrum Průhonice o.p.s.
- GEMINI eye clinic a.s., Zlín, Průhonice
- D. C. M. clinic, Hradec Králové
- FN Brno Transfusion & Tissue Department, Třebíč

• ...

Literature and references

- [1] Jamie M. Zeitzer, Derk-Jan Dijk, Richard E. Kronauer, Emery N. Brown and Charles A. Czeisler. Sensitivity of the human circadian pacemaker to nocturnal light: melatonin phase resetting and suppression. *Journal of Physiology.* 2000, 526.3, pp. 695–702.
- [2] Gooley, J. J. et al.: Spectral Responses of the Human Circadian System Depend on the Irradiance and Duration of Exposure to Light. Science, *Translational Medicine*. 2010, issue 31, ISSN 1946-6242, DOI 10.1126/scitranslmed. 3000741.
- [3] Dennis M. Dacey et al., Melanopsin-expressing ganglion cells in primate retina signal colour and irradiance and project to the LGN. *Letters to nature*. Nature vol. 433, 17 Feb. 2005. p. 749–754.
- [4] Thomas A. Wehr. Effect of Seasonal Changes in Daylength on Human Neuroendocrine Function. *Horm Res* 1998; 49. p. 118–124.
- [5] Terman MR, Wirz-Justice A. Chronotherapeutics for Affective Disorders: A Clinician's Manual for Light and Wake Therapy. 2nd edition, Kager 2013.
 ISBN 978-3-318-02090-8.

Thank you...

Antonín Fuksa

Further references

•www.ChBFT·CZ

Information in Czech, free booklet for download (English version available).

• fuksa@nasli.net